

**Materials List for:**

# Extraction of the EPP Component from the Surface EMG

Toshifumi Kumai<sup>1</sup>

<sup>1</sup>Graduate School of Oral Medicine, Matsumoto Dental University

URL: <https://www.jove.com/video/1653>

DOI: [doi:10.3791/1653](https://doi.org/10.3791/1653)

## Materials

Name	Company	Catalog Number	Comments
Analogue filter	Nihon Kohden	MEG6100	pre-,main-amplifier
A/D converter	National Instruments	PCI-MIO-16E-4	PCI board for computer
Connection interface	National Instruments	BNC-2090	8ch BNC adaptor
Disc electrodes	Nihon Kohden	NS-11	Ag/AgCl (Φ8mm)
Electrode past	Sanshin	SA-5	Semi-fluid carbon past
LabVIEW (Digital filter)	National Instruments	V. 8.5	Programming language
Chewing gum	Lotte	Green gum	Test food in clenching