

Materials List for:

Operant Sensation Seeking in the Mouse

Christopher M. Olsen¹, Danny G. Winder¹



¹Department of Molecular Physiology and Biophysics, Center for Molecular Neuroscience, Kennedy Center for Human Development, Vanderbilt University Medical Center

Correspondence to: Christopher M. Olsen at chris.olsen@vanderbilt.edu

URL: <https://www.jove.com/video/2292>

DOI: [doi:10.3791/2292](https://doi.org/10.3791/2292)

Materials

Name	Company	Catalog Number	Comments
Drug self-administration test package for mouse: extra-wide chamber and retractable levers	 Med Associates, Inc.	MED-307W-CT-D1	Levers are ultra-sensitive (require ~2 grams force) and are mounted 2.2 cm above the floor. Yellow stimulus lamps are mounted 2 cm above each lever.
Interface and software package	 Med Associates, Inc.	MED-SYST-16	This is the package for up to 16 chambers.