

Materials List for:

Ex Vivo Culture of Primary Human Fallopian Tube Epithelial Cells

Susan Fotheringham^{1,2}, Keren Levanon^{1,2,3}, Ronny Drapkin^{1,2,4}

¹Department of Medical Oncology, Dana-Farber Cancer Institute

²Harvard Medical School, Boston, MA

³Sheba Cancer Research Center, Chaim Sheba Medical Center

⁴Department of Pathology, Brigham and Women's Hospital

Correspondence to: Ronny Drapkin at ronny_drapkin@dfci.harvard.edu

URL: <https://www.jove.com/video/2728>

DOI: [doi:10.3791/2728](https://doi.org/10.3791/2728)

Materials

Name	Company	Catalog Number	Comments
Collagen from human placenta (Bornstein and Traub Type IV)	Sigma-Aldrich	C7521	
DMEM:F12 media	Cellgro	15-090-CV	
Ultrosor G	Crescent Chemical Company	67042	Add 20 ml of DMEM:F12 to stock bottle and use 10mls per 500 ml of media (2% final concentration)
Pen Strep	Invitrogen	15140-122	
BD Primaria culture plates	BD Biosciences	353802	
24 well Transwell Permeable supports	Corning	3470	Clear, 6.5 mm insert, 0.4 µm pore size, treated polyester
MEM (Minimal Essential Media)	Cellgro	10-010-CV	
Pronase	Roche Group	11459643001	
DNase	Sigma-Aldrich	DN25	
Sall2 antibody	Harvard - T. Benjamin Lab	Gift	1:20 dilution
Pax8 antibody	Proteintech	10336-1-AP	1: 1000 dilution