

**Materials List for:**

# A Noninvasive Hair Sampling Technique to Obtain High Quality DNA from Elusive Small Mammals

Philippe Henry<sup>1</sup>, Alison Henry<sup>1</sup>, Michael A. Russello<sup>1</sup>

<sup>1</sup>Department of Biology and Centre for Species at Risk and Habitat Studies,, University of British Columbia, Okanagan Campus

Correspondence to: Philippe Henry at [phenry@interchange.ubc.ca](mailto:phenry@interchange.ubc.ca)

URL: <https://www.jove.com/video/2791>

DOI: [doi:10.3791/2791](https://doi.org/10.3791/2791)

## Materials

### Hair snare:

- Rolls of clear packing tape (3M, St. Paul, MN, USA)
- Fishing line (10lb)
- Cryogenic tube (2ml)
- Sterile forceps
- Liquid Nitrogen Dewar (optional)

### DNA extraction:

- DNA IQ System (Cat. # DC6701; Promega) with the tissue and hair extraction Kit (Cat. # DC6740; Promega).
- MagneSphere Technology Magnetic Separation Stand (Cat. # Z5342; Promega)
- Micro centrifuge
- Hybridization oven or water bath
- 1.5 ml Eppendorf Tubes