

Materials List for:

Using an α -Bungarotoxin Binding Site Tag to Study GABA A Receptor Membrane Localization and Trafficking

Megan L. Brady¹, Charles E. Moon¹, Tija C. Jacob¹

¹Pharmacology & Chemical Biology Department, University of Pittsburgh School of Medicine

Correspondence to: Tija C. Jacob at tcj11@pitt.edu

URL: <https://www.jove.com/video/51365>

DOI: [doi:10.3791/51365](https://doi.org/10.3791/51365)

Materials

Name	Company	Catalog Number	Comments
P3 Primary Cell 4D kit	Lonza	V4XP-3024	
α -Bungarotoxin, Alexa Fluor 647 conjugate	Molecular Probes	B35450	
α -Bungarotoxin, Alexa Fluor 488 conjugate	Molecular Probes	B13422	
α -Bungarotoxin, Alexa Fluor 594 conjugate	Molecular Probes	B13423	
(+)-Tubocurarine chloride	Tocris	2820	
Mounting medium	DAKO	cS703	
DPBS no calcium, magnesium	Invitrogen	14190-136	
Poly-D-lysine	Sigma	P6407	
Gephyrin antibody	Synaptic Systems	147 011	1:300
VIAAT/VGAT antibody	Synaptic Systems	131 002	1:1,000
anti GFP	Life Technologies	A6455	1:2,000
Glass bottom tissue culture dish	MatTek Corporation	P35G-1.5-14-C - Mattek Dishes	
Coverslips (round cover glass), #1 thickness, 12 mm	Warner Instruments	640-0702	