

Materials List for:

Zebrafish Keratocyte Explants to Study Collective Cell Migration and Reepithelialization in Cutaneous Wound Healing

Jose L. Rapanan¹, Agnes S. Pascual¹, Chandana K. Uppalapati², Kimbal E. Cooper¹, Kathryn J. Leyva², Elizabeth E. Hull¹

¹Biomedical Sciences Program, Midwestern University

²Department of Microbiology & Immunology, Midwestern University

Correspondence to: Elizabeth E. Hull at ehullx@midwestern.edu

URL: <https://www.jove.com/video/52489>

DOI: [doi:10.3791/52489](https://doi.org/10.3791/52489)

Materials

Name	Company	Catalog Number	Comments
35 mm Glass Bottom Culture Dishes (Poly-D Lysine Coated)	MatTek	P35GC-1.5-14-C	1.5 mm thickness is best for coverslip removal. 14 mm diameter provides more culture area.
Swiss Jewelers Style Forceps, Integra, Miltex 17-303	VWR	21909-460	We find that narrow, fine forceps work the best for scale removal.