

Materials List for:

Basic Methods for the Study of Reproductive Ecology of Fish in Aquaria

Kazuya Fukuda¹, Tomoki Sunobe¹

¹Laboratory of Fish Behavioral Ecology, Tateyama Station, Field Science Center, Tokyo University of Marine Science and Technology

Correspondence to: Kazuya Fukuda at kazuya.f47@gmail.com

URL: <https://www.jove.com/video/55964>

DOI: [doi:10.3791/55964](https://doi.org/10.3791/55964)

Materials

Name	Company	Catalog Number	Comments
Hand net	Nisso	AQ-17	Select for the target species size.
Polyethylene bag	San-U Fish Farm	8194	
Rubber band	ESCO Co. LTD.	78-0420-64	ø 80 mm x 6 mm
Oxygen cylinder	N/A	N/A	Oxygen cylinder for diving equipment suits.
Elbagin	Japan Pet Design Co. Ltd.	75950	Pafurazine F (provided from same company) is equivalent drug to Elbagin.
Polystyrene foam box	N/A	N/A	
Pipette	AS ONE	1-8625-04	
Rope	Mizukami Kinzoku Co. LTD.	95301601	
Weight	N/A	N/A	Weight for diving equipment suits.
Water tank	N/A	N/A	
Air pump	KOTOBUKI	4972814 062115	
Air stone	KOTOBUKI	4972814 232204	
2-Methylquinoline	Wako	170-00376	
Ethanol (99.5)	Wako	057-00456	
Visible implant elastomer tag kit	Northwest Marine Technology	N/A	http://www.nmt.us/products/vie/vie.shtml
Soft sponge	N/A	N/A	
PVC board	N/A	N/A	0.3 mm thickness is easy to use.
Petri dish	N/A	N/A	A large one, such as ø 160 mm and 30 mm depth, is convenient for the injection of the VIE tag.
Transparent acrylic board	N/A	N/A	
UVA filtered light	N/A	N/A	
PVC pipe	N/A	N/A	ø 5 cm
Waterproof sheet	SOMAR Corp.	3EKW03	The film for the plain copier.
Sand	N/A	N/A	
Stereo microscope	N/A	N/A	
Camera	N/A	N/A	