

Materials List for:

The Unpredictable Chronic Mild Stress Protocol for Inducing Anhedonia in Mice

Or Burstein¹, Ravid Doron^{1,2}

¹School of Behavioral Science, The Academic College Tel-Aviv-Yaffo

²Department of Education and Psychology, Open University

Correspondence to: Ravid Doron at raviddor@mta.ac.il

URL: <https://www.jove.com/video/58184>

DOI: [doi:10.3791/58184](https://doi.org/10.3791/58184)

Materials

Name	Company	Catalog Number	Comments
Heating lamp	Ikea	AA-19025-3	
Heating pillow	Sachs	EF-188B	
Mice restrainer			
Portable electronic balance (*.** g)			
Standard rubber stopper, size 5	Ancare	#5.5R	To avoid spillage during SPT
Straight open drinking tube (2.5")	Ancare	OT-100	To avoid spillage during SPT (insert drinking tube into rubber stopper)
2% sucrose solution			
50 mL conical centrifuge tube			For the SPT
Pre-adolescent (approximately 20-days old) ICR outbred mice	Envigo	Hsd:ICR (CD-1)	