

Materials List for:

In Vitro ELISA Test to Evaluate Rabies Vaccine Potency

Corinne Jallet¹, Noël Tordo^{1,2}

¹Unit Antiviral Strategies, OIE Reference Laboratories for Rift Valley Fever Virus & Crimean Congo Hemorrhagic Fever Virus, Institut Pasteur

²Institut Pasteur de Guinée, Conakry

Correspondence to: Noël Tordo at ntordo@pasteur.fr

URL: <https://www.jove.com/video/59641>

DOI: [doi:10.3791/59641](https://doi.org/10.3791/59641)

Materials

Name	Company	Catalog Number	Comments
Class II Biological Safety Cabinet	ThermoFisher Scientific	10445753	if titrating live virus
Clear Flat-Bottom Immuno Nonsterile 96-Well Plates, 400 µL, MAXISORP	ThermoFisher Scientific	439454	good for binding to the loaded antibody
Equip Labo Polypropylene Laboratory Fume Hood	ThermoFisher Scientific	12576606	for the preparation of sulfuric acid
Immunology Plate Strong Adsorption MAXISORP Flat Bottom Well F96	Dutscher	55303	good for binding to the loaded antibody
Microplate Sealing Tape(100 sheets)	ThermoFisher Scientific	15036	
Microplate single mode reader Sunrise	TECAN		
Microplate shaker-incubator	Dutscher	441504	
Microplate washer Wellwash	ThermoFisher Scientific	5165000	
Multichannel pipette (30-300 µL) 12 channels	ThermoFisher Scientific	4661180N	
Single Channel pipettes (Kit 2 : Finnpiettes F2 0.2-2 µL micro, 2-20 µL, 20-200 µL & 100-1000 µL)	ThermoFisher Scientific	4700880	