

Materials List for:

Making Patch-pipettes and Sharp Electrodes with a Programmable Puller

Austin L. Brown¹, Brandon E. Johnson², Miriam B. Goodman²

¹Department of Molecular and Cellular Physiology, Stanford University

²Department of Molecular and Cellular Physiology, Stanford University School of Medicine

Correspondence to: Miriam B. Goodman at mbgoodman@stanford.edu

URL: <https://www.jove.com/video/939>

DOI: [doi:10.3791/939](https://doi.org/10.3791/939)

Materials

Name	Type	Company	Catalog Number	Comments
Micropipette Puller	Instrument	Sutter Instrument Co.	P-97	Or similar instrument (e.g. Sutter P-87 or P-2000)
Glass Capillaries	Reagent	Sutter Instrument Co.	BF150-86-10	Or, similar capillary glass. To make filling the pipette easier, use a capillary with a glass filament.