

Materials List for:

Micro-dissection of Rat Brain for RNA or Protein Extraction from Specific Brain Region

Kin Chiu¹, Wui Man Lau¹, Ho Tak Lau¹, Kwok-Fai So¹, Raymond Chuen-Chung Chang¹

¹Laboratory of Neurodegenerative Diseases, Department of Anatomy, LKS Faculty of Medicine, The University of Hong Kong - HKU

Correspondence to: Raymond Chuen-Chung Chang at rccchang@hku.hk

URL: <https://www.jove.com/video/269>

DOI: [doi:10.3791/269](https://doi.org/10.3791/269)

Materials

Name	Company	Catalog Number	Comments
Microscope	Olympus Corporation		
Dissection tools			Metal plate, forceps, blade, etc. Individually wrapped with aluminum foil and heated at 180oC for eight hours to eliminate RNase. Made ice cold before use.