

Materials List for:

Targeting Olfactory Bulb Neurons Using Combined *In Vivo* Electroporation and Gal4-Based Enhancer Trap Zebrafish Lines

Kenric J. Hoegler¹, Martin Distel², Reinhard W. Köster³, John H. Horne¹

¹Department of Biology, Pace University

²Cellular and Molecular Medicine, University of California, San Diego

³Division of Cell Biology and Cell Physiology, Zoological Institute, Braunschweig University of Technology

Correspondence to: John H. Horne at jhorne@pace.edu

URL: <https://www.jove.com/video/2964>

DOI: [doi:10.3791/2964](https://doi.org/10.3791/2964)

Materials

Reagents:

- Agarose, low-melting-point (Sigma, A9414)
- Plasmid midi- or maxi-prep kit (e.g. Qiagen, #12643)
- Phenol red (Sigma, P3532)

Equipment:

- Capillary glass with filament; ID = 0.78 mm; OD = 1.0 mm (Sutter Instruments, #BF100-78-10)
- Grass SD-9 Stimulator (Grass Technology, Model SD-9)
- Grass platinum subdermal electrodes, straight needle (Grass Technology, #E2)
- Micro-loader pipette tips (Eppendorf, #X22703R)
- Pipette puller (e.g. Sutter, Model P-30)
- Pipette puller (e.g. Sutter, Model P-30)