

Materials List for:

Curtain Flow Column: Optimization of Efficiency and Sensitivity

Sercan Pravadali-Cekic¹, Danijela Kocic¹, Stanley Hua¹, Andrew Jones¹, Gary Dennis¹, Andrew Shalliker¹

¹School of Science and Health, University of Western Sydney

Correspondence to: Andrew Shalliker at r.shalliker@uws.edu.au

URL: <https://www.jove.com/video/53471>

DOI: [doi:10.3791/53471](https://doi.org/10.3791/53471)

Materials

Name	Company	Catalog Number	Comments
HPLC instrument			
Additional Pump			Required if 2 pump CF system set up is to be used.
Curtain Flow HPLC column	Thermo Fisher Scientific	Not Defined	Soon to be commercialized
Methanol	Any brand		HPLC Grade
PEEK tubing	Any brand		Various lengths and i.d.
PEEK tube cutter	Any brand		
Analytical Scale Balance	Any brand		
Stop watch	Any brand		
Eluent collection vessels	Any brand		1-2 ml sample vials can be used as eluent collection vessels
T-piece	Any brand		