

Materials List for

# Identification of RNAs Engaged in Direct RNA-RNA Interaction with a Long Non-Coding RNA

Audrey Jacq<sup>1</sup>, Denis Becquet<sup>1</sup>, Séverine Guillen<sup>1</sup>, Bénédicte Boyer<sup>1</sup>, Maria-Montserrat Bello-Goutierrez<sup>1</sup>, Jean-Louis Franc<sup>\*1</sup>, Anne-Marie François-Bellan<sup>\*1</sup>

<sup>1</sup>Faculté des Sciences Médicales et Paramédicales, Aix Marseille Univ, CNRS, INP

\* These authors contributed equally

## Corresponding Authors

**Jean-Louis Franc**

jean-louis.franc@univ-amu.fr

**Anne-Marie François-Bellan**

anne-marie.francois@univ-amu.fr

## Citation

Jacq, A., Becquet, D., Guillen, S., Boyer, B., Bello-Goutierrez, M.M., Franc, J.L., François-Bellan, A.M. Identification of RNAs Engaged in Direct RNA-RNA Interaction with a Long Non-Coding RNA. *J. Vis. Exp.* (173), e62475, doi:10.3791/62475 (2021).

## Date Published

July 9, 2021

## DOI

10.3791/62475

## URL

jove.com/video/62475

## Materials

Name	Company	Catalog Number	Comments
4'-Aminomethyltrioxsalen hydrochloride	Sigma	A4330	Crosslinker reagent
Bioruptor Plus	Diagenode	B01020001	Sonicator
Biotynilated probes	IDT		Oligonucleotide probes
CFX96 Real Time System	BioRad	4351107	qPCR apparatus
DNA Oligonucleotides	IDT		Primers for qPCR
Dynabeads My One	Thermo-Fisher	65001	Magnetic streptavidin beads
Formamide	Thermo-Fisher	15515-026	Formamide
iTaq Universal SYBR Green Supermix	BioRad	1725124	qPCR reagent
Proteinase K	Sigma	P2308	Proteinase K
RNA to DNA	Thermo-Fisher	4387405	Reverse transcription kit
RNA XS purification kit	Macherey-Nagel	740902	RNA purification kit
RNaseOUT	Thermo-Fisher	10777-019	RNase inhibitor
Tube Rotator	Stuart	SB2	Eppendorf tube rotator
UV Stratalinker 1800	Stratagene	#400072	UV crosslinker