

Materials List for

# Rat Model of Normothermic Ex-Situ Perfused Heterotopic Heart Transplantation

Mukhammad Kayumov<sup>\*1,2,3</sup>, In Seok Jeong<sup>\*2,3</sup>, Dowan Kim<sup>2,3</sup>, Yujin Kwak<sup>2</sup>, Francis O. Obiweluzor<sup>2,3</sup>, Namsik Yoon<sup>4</sup>, Hyo Sin Kim<sup>5</sup>, Hwa Jin Cho<sup>3,6</sup>

<sup>1</sup>Department of Medical Science, Chonnam National University Graduate School <sup>2</sup>Department of Thoracic and Cardiovascular Surgery, Chonnam National University Hospital and Medical School <sup>3</sup>Cardiovascular and Respiratory Research Team, Chonnam National University Hospital <sup>4</sup>Department of Internal Medicine, Chonnam National University Hospital and Medical School <sup>5</sup>Departments of Surgery, Chonnam National University Hospital and Medical School <sup>6</sup>Department of Pediatrics, Chonnam National University Children's Hospital, and Medical School

\* These authors contributed equally

## Corresponding Authors

**Hyo Sin Kim**  
chhj98@gmail.com

**Hwa Jin Cho**  
gideon0504@gmail.com

## Citation

Kayumov, M., Jeong, I.S., Kim, D., Kwak, Y., Obiweluzor, F.O., Yoon, N., Kim, H.S., Cho, H.J. Rat Model of Normothermic Ex-Situ Perfused Heterotopic Heart Transplantation. *J. Vis. Exp.* (194), e64954, doi:10.3791/64954 (2023).

## Date Published

April 21, 2023

## DOI

10.3791/64954

## URL

jove.com/video/64954

## Materials

Name	Company	Catalog Number	Comments
AES active evacuation system	Smiths medical	PC-6769-51A	Utilize CO2 and excess isoflurane
Anesthesia machine	Smiths medical	PC-8801-01A	Mixes isoflurane and oxyegn and delivers to animal
B20 patient monitor	GE medical systems	B20	to observe mean aortic pressure and temperature
Homeothermic Monitoring System	Harvard apparatus	55-7020	To monitor and maintain animal's temperature
Micro-1 Rat oxygenator	Dongguan Kewei medical instruments	Micro-MO	For gas exchange in the langendorff circuit
Micropuncture introducer Set	COOK medical	G48007	for delivering cardioplegic solution to the arch through the abdominal aorta
Microscope	Amscope	MU1403	For zooming surgical field (Recipient)
Surgical loupe	SurgiTel	L2S09	For zooming surgical field (Donor)
Syringe pump	AMP all	SP-8800	To deliver cardioplegic solution
Transonic flow sensor	Transonic	ME3PXL-M5	Perfusion circuit flow sensor
Transonic tubing flow module	Transonic	TS410	flow acquiring system
Watson - Marlow pumps	Harvard apparatus	010.6131.DAO	Peristaltic pump used for recirculate perfusate
WBC-1510A	JEIO TECH	E03056D	Heating bath
Sprague-Dawley rats	Samtako Bio Korea Co., Ltd., Osan City Korea		
<b>Medications</b>			
BioHAnce Gel Eye Drops	SENTRIX Animal care		wet ointments for eye
Cefazolin	JW pharmaceutical		For prophylaxis
Custodiol	DR, FRANZ KOHLER CHEMIE GMBH		For heart harvesting

---

Diclofenac	Myungmoon Pharm. Co. Ltd		For pain control
Heparin	JW pharmaceutical		Anticoagulant
Insulin	JW pharmaceutical		hormon therapy
Saline	JW pharmaceutical		For hydration therapy