

Materials List for:

Chromosomal Spread Preparation of Human Embryonic Stem Cells for Karyotyping

Priscila B. Campos¹, Rafaela C. Sartore¹, Stacie N. Abdalla¹, Stevens K. Rehen¹

¹Institute of Biomedical Sciences, Federal University of Rio De Janeiro-UFRJ

Correspondence to: Stevens K. Rehen at srehen@anato.ufrj.br

URL: <https://www.jove.com/video/1512>

DOI: [doi:10.3791/1512](https://doi.org/10.3791/1512)

Materials

Name	Company	Catalog Number	Comments
15 ml centrifuge tube	Techno Plastic Products	91015	
Colcemid Karyo MAX	GIBCO, by Life Technologies	15212-012	
EDTA	Isofar	721	
Glacial acetic acid	Isofar	100	
Methanol	Isofar	208	
Potassium Chloride (KCl)	Merck & Co., Inc.	104.931.000	
Slide	Bioslide	7105-1	
Trypsin	Sigma-Aldrich	T4799	