

Materials List for

siRNA Encapsulation into Exosomes Using Electroporation: An Electric Pulse Based Technique to Load siRNA into Exosomes

1
?

URL: <https://www.jove.com/video/20566>
DOI: [doi:10.3791/20566](https://doi.org/10.3791/20566)

Materials

Name	Company	Catalog Number	Comments
Atto655-siRNA	Eurogentec	SQ-SIRNA	(Labelled-S) UGC-GCU-ACGAUC-GAC-GAU-G55; (UnlabelledAS) CAU-CGU-CGA-UCG-UAGCGC-A55.
Glass pipettes	Fisher Scientific	1156-6963	
Amaxa Nucleofector I	Lonza	Nucleofector I	With Amaxa's Nucleofector kits
Citric acid buffer with EDTA			Mix 0.1954 g citric acid and 0.2087 g disodium phosphate in 50 mL of deionised water. Add EDTA to 0.1 mM. Adjust pH to 4.4.
Microfuge tubes	Starlab	S1615-5500	1.5 ml