

Materials List for:

# A Murine Model of Group B *Streptococcus* Vaginal Colonization

Kathryn A. Patras<sup>1</sup>, Kelly S. Doran<sup>1,2</sup>

<sup>1</sup>Department of Pediatrics, Division of Host-Microbe Systems & Therapeutics, University of California San Diego School of Medicine

<sup>2</sup>Department of Biology and Center for Microbial Sciences, San Diego State University

Correspondence to: Kelly S. Doran at [kdoran@mail.sdsu.edu](mailto:kdoran@mail.sdsu.edu)

URL: <https://www.jove.com/video/54708>

DOI: [doi:10.3791/54708](https://doi.org/10.3791/54708)

## Materials

Name	Company	Catalog Number	Comments
Sesame oil	Sigma Aldrich	S3547-250ML	
β-Estradiol	Sigma Aldrich	E8875-1G	CAUTION: Wear appropriate PPE. β-estradiol can be absorbed through the skin and mucosal surfaces.
200 µl gel loading pipette tips	USA Scientific	1252-0610	
Urethro-genital, sterile, calcium alginate swabs	Puritan	25-801 A 50	
CHROMagar StrepB	DRG International	SB282	
Todd Hewitt Broth	Hardy Diagnostics	7161C	
18 G, 1.5 inch needles	BD	305199	
26 G, 0.5 inch needles	BD	305111	
10 ml syringes	BD	309604	
1 ml syringes	BD	309659	
0.45 µm PVDF syringe filters	Whatman	6900-2504	
Dulbecco's Phosphate-Buffered Salt Solution 1x	Corning	21-031-CV	