

Materials List for:

Application of Consistent Massage-Like Perturbations on Mouse Calves and Monitoring the Resulting Intramuscular Pressure Changes

Naoyoshi Sakitani^{1*}, Takahiro Maekawa^{1*}, Kumiko Saitou^{1,2}, Katsuhiko Suzuki³, Shuhei Murase^{1,4}, Masakuni Tokunaga¹, Daisuke Yoshino⁵, Keisuke Sawada⁶, Atsushi Takashima⁷, Motoshi Nagao¹, Toru Ogata¹, Yasuhiro Sawada^{1,8}

¹Department of Rehabilitation for Motor Functions, National Rehabilitation Center for Persons with Disabilities

²Graduate School of Sport Sciences, Waseda University

³Faculty of Sport Sciences, Waseda University

⁴Department of Orthopaedic Surgery, Graduate School of Medicine, The University of Tokyo

⁵Frontier Research Institute for Interdisciplinary Sciences, Tohoku University

⁶University of Cincinnati College of Medicine

⁷Department of Assistive Technology, National Rehabilitation Center for Persons with Disabilities

⁸Department of Clinical Research, National Rehabilitation Center for Persons with Disabilities

*These authors contributed equally

Correspondence to: Yasuhiro Sawada at ys454-ind@umin.ac.jp

URL: <https://www.jove.com/video/59475>

DOI: [doi:10.3791/59475](https://doi.org/10.3791/59475)

Materials

Name	Company	Catalog Number	Comments
Aluminum wire	DAISO JAPAN	B028	An aluminum wire is used to avoid escaping restriction by the wire
Blood pressure telemeter	Millar	SPR-671	A blood pressure telemeter is used to measure intramuscular pressure.
DAPI	Thermo Fisher Scientific	D1306	DAPI is a fluorescent probe which is commonly used to stain DNA for fluorescent microscopy.
Goat anti-rabbit Alexa Fluor 488 (Dilution ratio, 1:500)	Invitrogen	A11034	Antibody for immunohistochemical staining.
Goat anti-rat Alexa Fluor 568 (Dilution ratio, 1:500)	Invitrogen	A11077	Antibody for immunohistochemical staining.
ImageJ	NIH	N/A	Analysis software for image
LabChart8	ADInstruments		Analysis software for acquiring biological signals.
Prolong gold	Thermo Fisher Scientific	P36930	Prolong gold is for mounting stained samples.
Protein Block Serum-Free	Dako	X090930-2	For blocking non-specific background staining in immunohistochemical procedures.
Rat monoclonal anti-laminin-2 antibody (Dilution ratio, 1:1000)	Sigma Aldrich	L0663	Antibody for immunohistochemical staining.
Rat monoclonal anti-F4/80 antibody (Dilution ratio, 1:500)	Abcam	ab6640	Antibody for immunohistochemical staining.
Rabbit polyclonal anti-MCP-1 antibody (Dilution ratio, 1:1000)	Abcam	ab25124	Antibody for immunohistochemical staining.
Rabbit polyclonal anti-TNF- α antibody (Dilution ratio, 1:1000)	Abcam	ab66579	Antibody for immunohistochemical staining.
Surgical tape	3M Japan	1530EP-0	Surgical tape is used to restrict joint movement.