

Materials List for:

Infecting Mice with *Malassezia* spp. to Study the Fungus-Host Interaction

Florian Sparber¹, Salomé LeibundGut-Landmann¹

¹Section of Immunology, Vetsuisse Faculty, University of Zürich

Correspondence to: Florian Sparber at florian.sparber@uzh.ch, Salomé LeibundGut-Landmann at salome.leibundgut-landmann@uzh.ch

URL: <https://www.jove.com/video/60175>

DOI: [doi:10.3791/60175](https://doi.org/10.3791/60175)

Materials

Name	Company	Catalog Number	Comments
Agar	Sigma-Aldrich	A1296-1KG	
Attane Isoflurane	Piramal Healthcare	-	
Biosafety cabinet (BSC) FASTER Ultra Safe	DASIT GROUP	TEC 5594	BSL2 certified
Centrifuge	Eppendorf	5415D	compatible with 2ml Eppendorf tubes
Dessicated Ox-bile	Sigma-Aldrich	70168-100G	
Eppendorf Tubes (2 ml)	Eppendorf	0030 120.094	
Glucose	Sigma-Aldrich	49159-5KG	
Glycerol (99 %)	Honeywell	10314830	
Heating pad	Eickenmeyer	648048	
Incubator Heraeus B20	Heraeus	412047753	BSL2 certified
Ketasol (100 mg)	Graeb AG	6680416	
Magnetic heating plate MR Hei-Standard	Heidolph Instruments	442-1355	
<i>Malassezia</i> spp.	ATCC	14522, 14521, 42132	
Malt extract	Sigma-Aldrich	70167-500G	
Multiply Biosphere Tubes (200 µl)	Sarstedt AG	7084211	Safelock
Native olive oil	-	-	commerc. available
Nonidet P40	Axon Lab	A1694,0250	
Oditest measurement devise	Kroepelin	S0247	range 0-5 mm
Oleic Acid	Sigma-Aldrich	75090-5ML	
Peptone	Oxoid	LP0037	
Petri dishes	Sarstedt AG	82.1473	
Phosphat buffered salt solution (PBS, 1x)	Amimed/Bioconcept	3-05F39	
Rompun (2 %)	Bayer	KP0BFHR	
Shaking incubator Infors Minitron	Infors	-	BSL2 certified
Spectrometer	Jenway	20308	optical density measurement at 600nm
Spectrometer Cuvettes	Greiner Bio-One	613101	
Stainless Steel balls (5mm)	ABF	KU.5G80 1.3541	
Syringes 1 ml Sub-Q	BD Bioscience	305501	
Tissue Lyzer II	Quiagen	85300	
Transpore Hypoallergic Tape	3M	1527-1	
Tween 40	Sigma-Aldrich	P1504-100ML	

Vitamin A Retinoli Palmitas Eye Cream	BAUSCH & LOMB		commerc. available
---------------------------------------	---------------	--	--------------------