

Materials List for:

Retrieval of Mouse Oocytes

Amanda R. Duselis¹, Paul B. Vrana¹

¹Dept. of Biological Chemistry, University of California, Irvine (UCI)

Correspondence to: Amanda R. Duselis at aduselis@uci.edu, Paul B. Vrana at pvrana@uci.edu

URL: <https://www.jove.com/video/185>

DOI: [doi:10.3791/185](https://doi.org/10.3791/185)

Materials

| Name | Type | Company | Catalog Number | Comments |
|--------------------------------------|---------|---------------------------------|----------------|---|
| Flushing needle | Tool | Zephytronics | ZT5-130-L | 30 gauge blunt |
| Forceps (2) | Tool | Roboz Surgical Instruments Co. | RS-5060 | |
| Peromyscus | Animal | Peromyscus Genetic Stock Center | | |
| Scissors | Tool | Roboz Surgical Instruments Co. | RS-5880 | |
| FHM HEPES Buffered Media | Reagent | Specialty Media | MR-025-d | |
| KSOM Embryo Culture Medium (5x10 ml) | Reagent | Specialty Media | MR-020P-5F | |
| Mouth Pipette | Tool | | | Pulled from glass capillaries and attached to rubber tubing |
| Oocyte holding dish | Tool | | | Glass container with concave divets |