

Materials List for:

Colon Ascendens Stent Peritonitis (CASP) - a Standardized Model for Polymicrobial Abdominal Sepsis

Tobias Traeger¹, Pia Koerner¹, Wolfram Kessler¹, Katharina Cziupka¹, Stephan Diedrich¹, Alexandra Busemann¹, Claus-Dieter Heidecke¹, Stefan Maier¹

¹Department of Surgery, University of Greifswald

*These authors contributed equally

Correspondence to: Tobias Traeger at traeger@uni-greifswald.de

URL: <https://www.jove.com/video/2299>

DOI: [doi:10.3791/2299](https://doi.org/10.3791/2299)

Materials

Name	Company	Catalog Number	Comments
Narcotic fluid			Mix 0.5 mL Rompun with 4.0 mL Ketanest, and add 5.5 mL of saline solution (0.9%).For complete anesthesia, inject intraperitoneally 10 µL/g body weight
Rompun (20 mg/ml Xylazine)	Bayer AG		
Ketanest (25mg/ml Esketamine)	Pfizer Pharma GmbH	647028001E	
BD Venflon 18GA (1,2x45mm, 80ml/min)	BD Biosciences	391457	
BD Venflon 16GA (1,7x45mm, 180ml/min)	BD Biosciences	391455	
BD Venflon 14GA (2,0x45mm, 270ml/min)	BD Biosciences	391456	
Catgut Polyester 4/0 white, non absorbable	Catgut GmbH	17218113	
Mariderm black 7/0, non absorbable	Catgut GmbH	18104900	