

Materials List for:

Prediction of HIV-1 Coreceptor Usage (Tropism) by Sequence Analysis using a Genotypic Approach

Saleta Sierra¹, Rolf Kaiser¹, Nadine Lübke¹, Alexander Thielen², Eugen Schuelter¹, Eva Heger¹, Martin Däumer³, Stefan Reuter⁴, Stefan Esser⁵, Gerd Fätkenheuer⁶, Herbert Pfister¹, Mark Oette⁷, Thomas Lengauer²

¹Institute of Virology, University of Cologne

²Max Planck Institute for Informatics

³Institute for Immune genetics

⁴Department of Gastroenterology, Hepatology and Infectiology, University of Duesseldorf

⁵Department of Dermatology, University of Essen

⁶Department of Internal Medicine, University of Cologne

⁷Augustinerinnen Hospital

Correspondence to: Rolf Kaiser at rolf.kaiser@uk-koeln.de

URL: <https://www.jove.com/video/3264>

DOI: [doi:10.3791/3264](https://doi.org/10.3791/3264)

Materials

Name	Company	Catalog Number	Comments
MagNA Pure Compact Nucleic Acid Isolation Kit I	Roche Group	04 802 993 001	
MagNA Pure Compact System	Roche Group	03731146001	
T3000 Thermocycler	Biometra	050-723	
One Step RT-PCR Kit	Qiagen	210215	
HotStarTaq Master Mix Kit	Qiagen	203445	
QIAquick PCR Purification Kit	Qiagen	28106	
Exonuclease I (Exo I)	Fermentas	EN0582	
FastAP Thermosensitive Alkaline Phosphatase	Fermentas	EF0651	
BigDye Terminator v3.1 Cycle Sequencing Kit	Applied Biosystems	4337457	
Sephadex G-50 Superfine	GE Healthcare	17-0041-01	
ABI Prism 3130 XL capillary sequencer	Applied Biosystems	3130XL	