

**Materials List for:**

**Habituation and Prepulse Inhibition of Acoustic Startle in Rodents**

Bridget Valsamis<sup>1</sup>, Susanne Schmid<sup>1</sup>




<sup>1</sup>Anatomy and Cell Biology, Schulich School of Medicine and Dentistry, University of Western Ontario

Correspondence to: Susanne Schmid at [Susanne.schmid@schulich.uwo.ca](mailto:Susanne.schmid@schulich.uwo.ca)

URL: <https://www.jove.com/video/3446>

DOI: [doi:10.3791/3446](https://doi.org/10.3791/3446)

**Materials**

Name	Company	Catalog Number	Comments
Startle box package	 Med Associates, Inc.	MED-ASR-PRO1-ADD ( <a href="http://www.med-associates.com/startle/startle.htm">http://www.med-associates.com/startle/startle.htm</a> )	Includes hardware & software
Animal Holder	 Med Associates, Inc.	ENV-264A ( <a href="http://www.med-associates.com/startle/startle.htm#animal">http://www.med-associates.com/startle/startle.htm#animal</a> )	Other sizes and types also available
USB Sound Pressure Level Measurement Package	 Med Associates, Inc.	ANL-929A-PC ( <a href="http://www.med-associates.com/behavior/audio/generator.htm#anl929a">http://www.med-associates.com/behavior/audio/generator.htm#anl929a</a> )	For calibration