

Materials List for:

Determination of Mammalian Cell Counts, Cell Size and Cell Health Using the Moxi Z Mini Automated Cell Counter

Gregory M. Dittami¹, Manju Sethi¹, Richard D. Rabbitt², H. Edward Ayliffe¹

¹Orflo Technologies



²University of Utah

Correspondence to: Gregory M. Dittami at greg@orflo.com

URL: <https://www.jove.com/video/3842>

DOI: [doi:10.3791/3842](https://doi.org/10.3791/3842)

Materials

Name	Company	Catalog Number	Comments
ORFLO Diluent	 ORFLO	MXA006	
Cassette Dispenser	ORFLO		Stores up to 25 cassettes for convenient dispensing
USB Cable	ORFLO		Connects instrument to PC/Mac or power adapter
Power Adapter (US and EU models only)	ORFLO		Connects USB cable to an AC outlet
USB Flash Drive	ORFLO		Stores Moxi Z software and user manual
Calibration Check Beads	ORFLO		
Electronic Calibration Cassette	 ORFLO		Electronic cassette for verifying proper system operation and calibration