

Materials List for:

Dynamic Visual Tests to Identify and Quantify Visual Damage and Repair Following Demyelination in Optic Neuritis Patients

Noa Raz¹, Michal Hallak¹, Tamir Ben-Hur¹, Netta Levin¹

¹Department of Neurology, The Agnes Ginges Center for Human Neurogenetics, Hadassah Hebrew-University Medical Center

Correspondence to: Netta Levin at netta@hadassah.org.il

URL: <https://www.jove.com/video/51107>

DOI: [doi:10.3791/51107](https://doi.org/10.3791/51107)

Materials

Name	Company	Catalog Number	Comments
Personal computer, including laptops			The OFM software runs best on Mac or on Windows 7 (or higher) PC. The Stereo software runs on every personal computer.
Monitor specification			Size: at least 15inch, Color: at least 16 bit
The OFM and the Stereo softwares	These are self-developed softwares		Researchers & physicians who are interested in these softwares may contact us at: fmri-hadassah.org
Red/Cyan 3D glasses (We had tested the Stereo software on the two following 3D glasses):	Nvidia & American Paper Optics	3D Vision Ultimate Anaglyph 3D Glasses & Pro X Style Red/Cyan 3D Glasses for Movies and Games on Flat Screens	
Performance on our stereo task was compared to performance on the standard Randot stereo test	Stereo Optical Co.	Randot SO-002	