

Materials List for:

Treating Clinical Depression with Repetitive Deep Transcranial Magnetic Stimulation Using the Brainsway H1-coil

David Feifel¹, Katherine Pappas¹

¹Department of Psychiatry, University of California, San Diego

Correspondence to: David Feifel at dfeifel@mail.ucsd.edu

URL: <https://www.jove.com/video/53858>

DOI: [doi:10.3791/53858](https://doi.org/10.3791/53858)

Materials

Name	Company	Catalog Number	Comments
TMS Neurostimulator	Magstim	N/A	Magstim Rapid or Magstim Rapid ²
TMS H-coil	Brainsway	N/A	
TMS Cooling System	Brainsway	N/A	
TMS positioning System	Brainsway	N/A	
Mobile Medical Cart	Brainsway	N/A	The mobile medical cart contains the cooling system, the controller and the stimulator,
Personal cap and associated fixings	Brainsway	N/A	and enables mobility of the system inside the treatment area.
Ear plugs	Any	N/A	Rating of at least 30 dB noise reduction
Pillow	Any	N/A	Used for positioning arm when determining the motor threshold
dTMS instruction manual	Brainsway	N/A	