

Materials List for:

Visualization of Ambient Mass Spectrometry with the Use of Schlieren Photography

Gregory T. Winter¹, Joshua A. Wilhide¹, William R. LaCourse¹

¹Department of Chemistry and Biochemistry, University of Maryland, Baltimore County

Correspondence to: William R. LaCourse at lacourse@umbc.edu

URL: <https://www.jove.com/video/54195>

DOI: [doi:10.3791/54195](https://doi.org/10.3791/54195)

Materials

Name	Company	Catalog Number	Comments
Flashlight	EAGTAC	D25A Ti	or equivalent
Spherical Concave Mirror	Anchor Optics	27633	
Rebel EOS T2i	Canon	4462B001	or equivalent
300 mm telephoto lens	Canon	6473A003	or equivalent
Direct Sample Analysis (DSA) Ionization Source	PerkinElmer	MZ300560	or equivalent
Sq 300 MS with SQ Driver Software	PerkinElmer	N2910801	or equivalent
Ring Stand	Fisher Scientific	11-474-207	or equivalent
Laser Pointer	Apollo	MP1200	or equivalent
razor blade	Blue Hawk	34112	or equivalent
small drill bit #73	CML Supply	503-273	or equivalent
Protractor	Sterling	582	or equivalent
Hose Clamp	Trident	720-6000L	or equivalent