

Materials List for:

How to Use the H1 Deep Transcranial Magnetic Stimulation Coil for Conditions Other than Depression

Aron Tendler¹, Yiftach Roth¹, Noam Barnea-Ygael², Abraham Zangen²

¹Brainsway LTD

²Department of Life Sciences, Ben Gurion University of the Negev

Correspondence to: Aron Tendler at aronte@brainsway.com

URL: <https://www.jove.com/video/55100>

DOI: [doi:10.3791/55100](https://doi.org/10.3791/55100)

Materials

Name	Company	Catalog Number	Comments
dTMS System	Brainsway		Includes H1 coil, positioning arm, cart, stimulator, cooling system
Patient Caps	Brainsway		Includes blue caps with rulers
Ear plugs			Rated to 30 dB