

Materials List for:

Sample Preparation for Probe Electrospray Ionization Mass Spectrometry

Sen Takeda¹, Kentaro Yoshimura¹, Hiroshi Tanihata²

¹Department of Anatomy and Cell Biology, University of Yamanashi School of Medicine

²Shimadzu Corporation

Correspondence to: Sen Takeda at stakeda@yamanashi.ac.jp

URL: <https://www.jove.com/video/59942>

DOI: [doi:10.3791/59942](https://doi.org/10.3791/59942)

Materials

Name	Company	Catalog Number	Comments
5-Fluoro-2'-deoxyuridine (5-FdU)	Sigma-Aldrich	F8791-25MG	25mg
disposable biopsy punch (Trepan)	kai Europa GmbH	BP-30F	bore size 3mm
ethanol	nacalai tesque	14710-25	extra pure reagent
LabSolutions	Shimadzu		ver. 5.96, Data analyzer
micropestle	United Scientific Supplies	S13091	
microtube	Treff	982855	0.5 mL clear
PESI-MS (Direct Probe Ionization-MS)	Shimadzu	DPiMS-2020	Mass spectrometer equipped with PESI
PPGT solution	Shimadzu	ND	Attached to DPiMS-2020