

Materials List for

## Cigarette Smoke Exposure in Mice using a Whole-Body Inhalation System

Daniel E. Morales-

Mantilla<sup>1,2</sup>, Xinyan Huang<sup>3,7</sup>, Philip Erice<sup>1,3</sup>, Paul Porter<sup>4</sup>, Yun Zhang<sup>1,5</sup>, Mary Figueroa<sup>6</sup>, Joya Chandra<sup>6</sup>, Katherine Y. King<sup>2</sup>, Farrah Kheradmand<sup>4,8</sup>, Antony Rodríguez<sup>3,8</sup>

<sup>1</sup>Program in Immunology, Baylor College of Medicine <sup>2</sup>Department of Pediatrics, Section of Infectious Diseases, Baylor College of Medicine <sup>3</sup>Department of Medicine - Immunology Allergy and Rheumatology, Baylor College of Medicine <sup>4</sup>Department of Medicine, Pulmonary, Critical Care, Sleep Medicine, Baylor College of Medicine <sup>5</sup>Department of Pathology and Immunology, Baylor College of Medicine <sup>6</sup>Department of Pediatrics, Research and Department of Epigenetics and Molecular Carcinogenesis, The University of Texas MD Anderson Cancer Center <sup>7</sup>Division of Pulmonary and Critical Care Medicine, The First Affiliated Hospital of Sun Yat-sen University <sup>8</sup>Center for Translational Research on Inflammatory Diseases (CTRID), Michael E. DeBakey VA Medical Center

Corresponding Author	Citation	Citation		
Antony Rodríguez	Morales-Mantilla, D.E., Huang,	Morales-Mantilla, D.E., Huang, X., Erice, P., Porter, P., Zhang, Y., Figueroa, M.,		
antonyr@bcm.edu	Chandra, J., King, K.Y., Kheradmand, F., Rodríguez, A. Cigarette Smoke Exposure in Mice using a Whole-Body Inhalation System. <i>J. Vis. Exp.</i> (), e61793, doi:10.3791/61793 (2020).			
Date Published	DOI	URL		
October 22, 2020	10.3791/61793	jove.com/video/61793		

## **Materials**

Name	Company	Catalog Number	Comments
1 in fastener	Lowes	756990	
1/4 in Barbed Y connector	VWR	89093-282	
1/4 in straight tubing connector	VWR	62866-378	
1/8 hex nipple	Lowes	877221	
1/8 in threaded coupling fitting	Lowes	877208	
1/8 in threaded male adapter nipple fitting	Lowes	877243	
10/32 (M) threaded straight connector	Bimba	EB60	
3/4 in 90-degree elbow CPVC fitting	Lowes	22643	
3/4 in chlorinated polyvinyl chloride (CPVC) pipe	Lowes	23814	
3/4 in CPVC cap	Lowes	23773	
3/4 in CPVC Drip irrigation female adapter	Lowes	194629	
3/4 in diameter CPVC male adapter	Lowes	23766	
8.5 L airtight container with lid (11.25in x 7.75in x 6 in)	Komax	N/A	Listed as "Komax Biokips Large Bread Box   (280-oz) Large Storage Container"
Glass drain tube (1.75 in diameter x 8 in length)	KIMAX	6500	
Isonic Solenoid Valves	Bimba	V2A02-AW1	
Marlboro Red 100's	Marlboro	N/A	
Oxygen swivel barbed connector	Global Medical Solutions	RES002	



Panasonic Timer LT4H-W	Panasonic	LT4HW	Item was built-in the valve controller by Shepherd Controls & Associates
Pressure regulator	Allied Electronics and Automation	70600552	Also listed as "Norgren R07-100- RGKA"
Rubber stopper # 1 (one hole)	VWR	59581-163	
Rubber stopper # 8.5 (one hole)	VWR	59581-389	
Scireq inExpose system	Scireq and Emka Technologies	N/A	Commercial system used for comparison with our DIY WBIS
Straight barbed fitting (8mm opening)	VWR	10028-872	
Thread Sealant tape	Lowes	1184243	
Threaded port adaptor	Bimba	P1SA1	
Timeter Aridyne 2000 Medical Air Compressor	MFI Medical	AHC-TE20	
Timeter flowmeter	Allied Healthcare Products	15006-03YP2	Also listed as "Puritan Air Meter"
Valve Control system	Shepherd Controls and Associates	N/A	Company custom designed the valve control system for this model.
Vinyl pipes	Vitality Medical	RES3007	