

Materials List for:

C. elegans Gonad Extrusion: A Rapid Dissection Technique

1
?

URL: <https://www.jove.com/t/20131>

DOI: [doi:10.3791/20131](https://doi.org/10.3791/20131)

Materials

Name	Company	Catalog Number	Comments
Flat bottom glass watch dish	Agar Scientific	AGL4161	We use glass because dissected gonads often stick to plastic
25 G needles	BD PrecisionGlide	305122	
Syringes (could be 1 or 5 mL)	BD Syringes	1 mL: BD 309659	
Microscope slides (25 x 75 x 1.0 mm)	Fisherbrand	12-550-343	
Clinical bench top centrifuge			
*M9 solution			
Levamisole	Sigma	L9756	
*M9 Buffer Recipe			3 g KH ₂ PO ₄ , 6 g Na ₂ HPO ₄ , 5 g NaCl, 1 mL 1 M MgSO ₄ , H ₂ O to 1 L.