

Materials List for:

***In vivo* Visualization of Synaptic Vesicles Within *Drosophila* Larval Segmental Axons**

Michelle L. Kuznicki¹, Shermali Gunawardena¹

¹Department of Biological Sciences, SUNY-University at Buffalo

Correspondence to: Shermali Gunawardena at sg99@buffalo.edu

URL: <https://www.jove.com/video/2151>

DOI: [doi:10.3791/2151](https://doi.org/10.3791/2151)

Materials

Name	Company	Catalog Number	Comments
1 Pair Micro Dissection Scissors	Fine Science Tools		Spring scissors – 6mm blade. 0.125 mm tip diameter
2 Pair Forceps	Fisher Scientific		Fisherbrand Dissecting Micro-Adson serrated tip forceps
Box of Dissection Pins	Fine Science Tools		Minutien Pins 0.1mm diameter