

Materials List for:

Modulating Cognition Using Transcranial Direct Current Stimulation of the Cerebellum

Paul A. Pope¹

¹School of Psychology, University of Birmingham

Correspondence to: Paul A. Pope at p.pope@bham.ac.uk

URL: <https://www.jove.com/video/52302>

DOI: [doi:10.3791/52302](https://doi.org/10.3791/52302)

Materials

Name	Company	Catalog Number	Comments
Headset	Beyerdynamic	DT234 Pro	Ensure the microphone does not move from the participants mouth in between testing.
DC stimulator	Magstim	DC Stimulator Plus	Electrode placement is a critical success factor for tDCS efficacy
Stimulus presentation and recording software	www.neurobs.com	Presentation (Version 14.2)	Maintain participant-specific stimulus presentation rate between sessions in experiment one