

**Materials List for:**

**Ex Vivo Culture of Chick Cerebellar Slices and Spatially Targeted Electroporation of Granule Cell Precursors**

Michalina Hanzel<sup>1</sup>, Richard J.T. Wingate<sup>1</sup>, Thomas Butts<sup>2</sup>

<sup>1</sup>MRC Centre of Developmental Neurobiology, King's College London

<sup>2</sup>School of Biological and Chemical Sciences, Queen Mary, University of London

Correspondence to: Thomas Butts at [thomas.butts@qmul.ac.uk](mailto:thomas.butts@qmul.ac.uk)

URL: <https://www.jove.com/video/53421>

DOI: [doi:10.3791/53421](https://doi.org/10.3791/53421)

**Materials**

Name	Company	Catalog Number	Comments
Mcllwain tissue chopper	Mickle Laboratory Engineering Ltd		Cut at 300 µm for best results.
Basal Medium Eagle (Gibco)	Life Technologies	41010-026	
L-glutamine	Sigma	G7513	
penicillin/streptomycin	Sigma	P4333	
0.4 µm culture insert	Millipore	PICM0RG50	
TSS20 Ovodyne electroporator	Intracel	01-916-02	Use 3 x 10 V, 10 msec pulses for electroporation.