

Materials List for:

A Multi-hole Cryovial Eliminates Freezing Artifacts when Muscle Tissues are Directly Immersed in Liquid Nitrogen

Yizhong Huang¹, Maozhang He¹, Qingjie Zeng¹, Lin Li¹, Zhen Zhang¹, Junwu Ma¹, Yanyu Duan¹

¹State Key Laboratory for Pig Genetic Improvement and Production Technology, Jiangxi Agricultural University

Correspondence to: Yanyu Duan at yanyudian@hotmail.com

URL: <https://www.jove.com/video/55616>

DOI: [doi:10.3791/55616](https://doi.org/10.3791/55616)

Materials

Name	Company	Catalog Number	Comments
Cryostat Microtome	Leica	Leica CM1950	
Digital Microscope	Nikon	Nikon DS-U3	
Cryogenic Vial Plastic film			Designed by ourself
Liquid Nitrogen			Commonly-used
Scalpel			Commonly-used
10 cm forceps			Commonly-used
25 cm tweezer			Commonly-used
Safety glass			Commonly-used
Freezer gloves			Commonly-used