

Materials List for:

Left Atrial Stenosis Induced Pulmonary Venous Arterialization and Group 2 Pulmonary Hypertension in Rat

Ping Yu Xiong^{*1,2}, Shunsuke Baba^{*3,4}, Naritomo Nishioka^{3,5}, Yoshitaka Fujimoto³, Stephen L. Archer¹, Susumu Minamisawa^{1,4}

¹Department of Medicine, Queen's University

²Department of Biomedical and Molecular Sciences, Queen's University

³Department of Cell Physiology, The Jikei University School of Medicine

⁴Department of Pediatrics, The Jikei University School of Medicine

⁵Department of Cardiovascular Surgery, The Jikei University School of Medicine

*These authors contributed equally

Correspondence to: Susumu Minamisawa at sminamis@jikei.ac.jp

URL: <https://www.jove.com/video/58787>

DOI: [doi:10.3791/58787](https://doi.org/10.3791/58787)

Materials

Name	Company	Catalog Number	Comments
5-0 Prolene Suture	Johnson & Johnson - Ethicon	8725H	Polypropylene suture with HEMO-SEAL Technology
Anaesthesia Machine	Wakenyaku Co., Ltd.	BRTK-100A	Air pump and anaesthesia machine
Angiocatheter guidewire	Self-made		10 cm guidewire glued to a 1 cc syringe
Chest retractor	Natsume Seisakusho Co., Ltd.	F-2	
Chest tube 23G	Self-made		10 cc syringe attached to a 23G needle plus plastic tube
Curved forceps	Natsume Seisakusho Co., Ltd.	A-14	
Heating pad	Vivaria	MP-916-NV	Keep body temperature at 37 degree celsius
Horizon Ligating Clips	Teleflex	REF 003200	Size Medium-Large
Horizon Manual-Load Ligating Clip Applier For Medium-Large Size Horizon	Teleflex	REF 337085	Ligation Clips Angled Jaw, (20cm)
Needle holder	Natsume Seisakusho Co., Ltd.	MC-40	
Rodent Respirator	CWE Inc	SAR-830/P	Small animal ventilator
Scissors	Natsume Seisakusho Co., Ltd.	B-12	Straight scissors ideally with round tips
Straight forceps	Natsume Seisakusho Co., Ltd.	A-7	
Tongue depressor	Uchida Yoko Co., Ltd.	8-615-2417	Use the wide end