

**Materials List for:**

# Application of Voltage in Dynamic Light Scattering Particle Size Analysis

Tianyu Ren<sup>\*1</sup>, Emma J. Roberge<sup>\*2</sup>, John R. Csoros<sup>1</sup>, W. Rudolph Seitz<sup>1</sup>, Eva Rose M. Balog<sup>3</sup>, Jeffrey Mark Halpern<sup>2</sup>

<sup>1</sup>Department of Chemistry, University of New Hampshire

<sup>2</sup>Department of Chemical Engineering, University of New Hampshire

<sup>3</sup>Department of Chemistry and Physics, University of New England

\*These authors contributed equally

Correspondence to: Eva Rose M. Balog at [ebalog@une.edu](mailto:ebalog@une.edu), Jeffrey Mark Halpern at [jeffrey.halpern@unh.edu](mailto:jeffrey.halpern@unh.edu)

URL: <https://www.jove.com/video/60257>

DOI: [doi:10.3791/60257](https://doi.org/10.3791/60257)

## Materials

Name	Company	Catalog Number	Comments
N-Isopropylacrylamide	Tokyo Chemical Industry CO., LTD	I0401-500G	
1,4-Dioxane	Alfa Aesar	39118	
2,2"-Azobis(2-methylpropionitrile)	SIGMA-ALDRICH	441090-100G	
Cuvette	Malvern	DTS0012	
Dynamic Light Scattering	Malvern	Zetasizer NanoZS	
Ferrocenylmethyl methacrylate	ASTATECH	FD13136-1G	
Phthalimidomethyl butyl trithiocarbonate	SIGMA-ALDRICH	777072-1G	
Potentiostat	Gamry	Reference 600	