

## Materials List for

## **Assessing Functional Metrics of Skeletal Muscle Health in Human Skeletal Muscle Microtissues**

Heta Lad<sup>1,2</sup>, Brennen Musgrave<sup>1,2</sup>, Majid Ebrahimi<sup>1,2</sup>, Penney M. Gilbert<sup>1,2,3</sup>

<sup>1</sup>Institute of Biomedical Engineering, University of Toronto <sup>2</sup>Donnelly Centre for Cellular and Biomolecular Research, University of Toronto <sup>3</sup>Department of Cell and Systems Biology, University of Toronto

Corresponding Author	Citation			
Penney M. Gilbert	Lad, H., Musgrave, B., Ebrahim	Lad, H., Musgrave, B., Ebrahimi, M., Gilbert, P.M. Assessing Functional Metrics of		
Penney.Gilbert@utoronto.ca	Skeletal Muscle Health in Human Skeletal Muscle Microtissues. <i>J. Vis. Exp.</i> (), e62307, doi:10.3791/62307 (2021).			
Date Published	DOI	URL		
February 18, 2021	10.3791/62307	jove.com/video/62307		

## **Materials**

Name	Company	Catalog Number	Comments
0.9% Saline Solution, Sterile	House Brand	1010	10 mL aliquots of the solution are made and stored at 4°C
25G Needle	BD, Medstore, University of Toronto	2548-CABD305127	
6-Aminocaproic Acid, ≥99% (titration), Powder	Sigma - Aldrich	A2504-100G	A 50 mg / mL stock solution is generated by dissolving 5 mg of 6-aminocaproic acid powder in 100 mL of autoclaved, distilled water. The solution is vaccum filtered and 10 mL aliquots are stored at 4°C
6.35 mm ID Tubing	VWR	60985-528	
AB1167 Myoblast Cell Line	Institut de Myologie (Paris, France)		
Arbitrary Waveform Generator	Rigol	DG1022Z	
Basement Membrane Extract (Geltrex)	Thermo Fisher Scientific	A14132-02	Stored as aliquots of 50 µL or 100 µL at -80°C
Benchtop Vacuum Chamber	Sigma - Aldrich	D2672	
BNC to Aligator Clip Cable			Ordered from Amazon
Culture Plastics	Sarstedt		Includes culture plates, serological pipettes, etc
Dimethyl Sulfoxide	Sigma - Aldrich	D8418-250ML	
DPBS, Powder, No Calcium, No Magnesium	Thermo Fisher Scientific	21600069	
Dulbecco's Modified Eagle Medium (DMEM) (1X)	Gibco	11995-065	This is a high glucose DMEM with L-glutamine and sodium pyruvate
Fetal Bovine Serum	Fisher Scientific	10437028	
Fibrinogen from Bovine Plasma	Sigma - Aldrich	F8630-5G	Aliquots ranging from 7 - 10 mg of fibrinogen powder are made and stored at -20°C
Filtropur Syringe Filter, 0.22um Pore Size	Sarstedt	83.1826.001	
Horse Serum	Gibco	16050-122	
Human Recombinant Insulin	Sigma - Aldrich	91077C	Stock solution is 100X and made by dissolving 1 mg of human



			recombinant insulin in 1 mL of DMEM and 1 µL of NaOH 10N. Solution is filtered and stored as 1 mL aliquots at 4°C
Image Acquisition Software	Olympus	cellSens Dimension	
Image Processing Software	National Institutes of Health	ImageJ	
Isotemp Oven	Thermo Fisher Scientific	201	
Microscope	Olympus	IX83	
Microscope - Camera Mount	Labcam	Labcam for iPhone	Ordered from Amazon
Penicillin-Streptomycin (10,000 U/mL)	Gibco	15140-122	
Plastic Disposable Syringes, 1cc	BD	2606-309659	
Plastic Disposable Syringes, 50cc	BD	2612-309653	
Pluronic F-127, Powder, BioReagent	Sigma - Aldrich	P2443-250G	A 5% stock solution of pluronic acid is made by dissolving 5 g of pluronic acid powder in 100 mL of chilled, autoclaved, distilled water. The solution is vaccum filtered and 10 mL aliquots are stored at 4°C
Polydimethylsiloxane (Sylgard 184 Silicone Elastomer Kit)	Dow	4019862	Kits are also available at Thermo Fisher Scientific, Sigma - Aldrich, etc.
Polyurethane Negative Mold	In House		
Release Agent	Mann Release Technologies	200	
Rotary Vane Vacuum Pump	Edwards	A65401906	
Scalpel	Almedic, Medstore, University of Toronto	2586-M36-0100	
Single Edge Razor Blade	VWR	55411-050	
Skeletal Muscle Cell Basal Medium	Promocell	C-23260	30 mL aliquotes are generated and at stored at 4°C.
Skeletal Muscle Cell Growth Medium (Ready-to-use)	Promocell	C-23060	42 mL aliquots are generated and stored at 4°C.
Smartphone (iPhone)	Apple	SE	
Standard Duty Dry Vacuum Pump	Welch	2546B-01	
Sterilization Bag	Alliance	211-SCM2	
Thimble	Igege		Ordered from Amazon
Thrombin from human plasma	Sigma - Aldrich	T6884-250UN	100 units of thrombin is dissolved in 1 mL of a 0.1% BSA solution. 10 µL aliquots are prepared and stored at - 20°C.
Tin coated copper wire	Arco	B8871K48	Ordered from Amazon
Trypan Blue Solution, 0.4%	Thermo Scientific	15250061	
Trypsin-EDTA, 0.25%	Thermo FIsher Scientific	25200072	
Vacuum Chamber 2	SP Bel-Art	F42027-0000	