

Materials List for

# Robotic Left Hepatectomy using Indocyanine Green Fluorescence Imaging for an Intrahepatic Complex Biliary Cyst

Burak Görgec<sup>1,2</sup>, Babs M. Zonderhuis<sup>2,3</sup>, Marc G. Besselink<sup>1,2</sup>, Joris Erdmann<sup>1,2,3</sup>, Geert Kazemier<sup>2,3</sup>, Rutger-Jan Swijnenburg<sup>1,2,3</sup>

<sup>1</sup>Department of Surgery, Amsterdam UMC, location University of Amsterdam <sup>2</sup>Cancer Center Amsterdam <sup>3</sup>Department of Surgery, Amsterdam UMC, location Vrije Universiteit

## Corresponding Author

Rutger-Jan Swijnenburg

r.j.swijnenburg@amsterdamumc.nl

## Citation

Görgec, B., Zonderhuis, B.M., Besselink, M.G., Erdmann, J., Kazemier, G., Swijnenburg, R.J. Robotic Left Hepatectomy using Indocyanine Green Fluorescence Imaging for an Intrahepatic Complex Biliary Cyst. *J. Vis. Exp.* (184), e63265, doi:10.3791/63265 (2022).

## Date Published

June 24, 2022

## DOI

10.3791/63265

## URL

jove.com/video/63265

## Materials

Name	Company	Catalog Number	Comments
<i>Systems</i>			
Arietta V70 Ultrasound	Hitachi	-	The ultrasound system.
da Vinci Surgeon Console	IS	SS999	Used to control the surgical robot.
da Vinci Vision Cart	IS	VS999	The vision cart houses advanced vision and energy technologies and provides communications across da Vinci system components.
da Vinci Xi	IS	K131861	The surgical robot: 'patient side-cart'.
Robotic ultrasonography transducer	Hitachi	L43K	Used for intraoperative laparoscopic ultrasonography.
<i>Instruments</i>			
da Vinci Xi Endoscope with Camera, 8 mm, 30°	IS	470027	The camera of the da Vinci robot.
EndoWrist Fenestrated Bipolar Forceps	IS	470205	Used for dissection and coagulation.
EndoWrist HOT SHEARS	IS	470179	Used for cutting and coagulation.
EndoWrist Maryland Bipolar Forceps	IS	470172	Used for dissection.
EndoWrist Permanent Cautery Hook	IS	470183	Used for coagulation.
EndoWrist Medium-Large Clip Applier	IS	470327	Used for clipping with Weck Hem-o-lok medium-large polymer clip
EndoWrist Stapler 45 Instrument	IS	470298	Used for stapling
Vessel sealer	IS	480322	Used for vessel sealing and dividing.