

Materials List for:

Pavlovian Conditioned Approach Training in Rats

Christopher J. Fitzpatrick¹, Jonathan D. Morrow²

¹Neuroscience Graduate Program, University of Michigan

Correspondence to: Jonathan D. Morrow at jonmorro@med.umich.edu

URL: https://www.jove.com/video/53580 DOI: doi:10.3791/53580

Materials

| Name | Company | Catalog Number | Comments |
|--|----------------------|----------------|--|
| Standard Modular Operant Test Chamber with Modified Top for Rat | Med Associates, Inc. | ENV-008CT | Top has access slot for drug delivery leashes or other harness systems; comes with ENV-008CTC (Cover Insert) for access slot |
| Expanded PVC Sound Attenuating Cubicle | Med Associates, Inc. | ENV-022V | |
| PCI Operating Package for up to Sixteen Chambers | Med Associates, Inc. | MED-SYST-16 | Based on the number of operant chambers, MED-SYST-8 (PCI Operating Package for up to Eight Chambers) is also available |
| SmartCtrl | Med Associates, Inc. | DIG-716P2 | |
| Universal Cable, 25' (7.6 m) | Med Associates, Inc. | SG-210CB | Based on individual needs, SG-210CB-50 (Universal Cable, 50' [15.2 m]) is also available |
| 18" (45.7 cm) 3-Pin Mini-Molex | Med Associates, Inc. | SG-216A | Based on individual needs, SG-216A-2 (2' [61.0 cm] 3-Pin Mini-Molex) is also available |
| Power Cable, 25' (7.6 m) | Med Associates, Inc. | SG-210CP-25 | |
| Stainless Steel Grid Floor for Rat | Med Associates, Inc. | ENV-005 | Based on preference, ENV-005A (Stainless Steel Grid Floor for Mouse) is also available |
| Reusable Waste Pan | Med Associates, Inc. | ENV-007-P3 | |
| Filler Panel Package for Standard Modular Test Chamber | Med Associates, Inc. | ENV-008-FP | |
| Modular Pellet Dispenser for Rat, 45 mg | Med Associates, Inc. | ENV-203M-45 | |
| Pellet Receptacle, Trough Type | Med Associates, Inc. | ENV-200R2M | Based on preference, ENV-200R1M (Pellet Receptable, Cup Type) is also available |
| Retractable Lever | Med Associates, Inc. | ENV-112CM | Calibrated to 25 g; needs to be adjusted to 10 g before Pavlovian conditioned approach training |
| Retractable Lever Cue Light | Med Associates, Inc. | ENV-112CML | |
| House Light for Rat | Med Associates, Inc. | ENV-215M | |
| 100 mA Replacement Bulbs, Pack of 10 | Med Associates, Inc. | SG-800R | Red-colored; for use with ENV-215M (House Light for Rat) |
| Pellets, 45 mg, Primate Purified Diet, Banana Flavor | Bio-Serv | F0059 | 50,000/box |

²Department of Psychiatry, University of Michigan